


Wearables: Overview

JSM 2019

Wearables



source: North American Consumer Technographics Consumer Technology Survey, 2014

Wearables

Research



Consumer




What do wearables offer?

- Physical Activity
 - Steps or Activity Counts
 - Steps and Gait (temporal asymmetry, stride variability)
 - Energy Expenditure (calories, ...)
- Sleep
- Circadian Rhythmicity
- Electronic Diary/Ecological Momentary Assessment (1-2-4 per day)
- Heart Rate (ECG, bpm)
- Blood Glucose Monitoring
- Light, Temperature (Circadian markers)
- Voice (Mood, Progression of Disease)

Clinical applications

- Aging (BLSA, Health ABC, NHANES, UKBiobank, WHI)
- Dementia and AD (Sleep & Agitation)
- Cardiovascular: CHF, Afib, and post-surgery
- Multiple Sclerosis (Disability & Sleep)
- Mood Disorders
- Cancer: Fatigue and Sleep
- Diabetes (T2)
- Diabetes in babies (Nurture)
- Rehabilitation (METRC)


Accelerometers




- Detects acceleration in three orthogonal planes
- <https://www.youtube.com/watch?v=irjG9Y4NGnE>

Macro- and Micro-scale

- **Macro-scale** – summarized data (1 minute intervals)




- **Micro-scale** – raw accelerometry data collected (10Hz+)



Subsecond-level accelerometry data analysis

- **Vector Magnitude (VM):** 1-dimensional summary of 3-dimensional time-series


$$r(t) = \sqrt{x_1^2(t) + x_2^2(t) + x_3^2(t)}$$



Stage 1: Episode Detection


- Non-wear time
- Posture: sitting, lying, standing, driving, stairs climbing, ...
- Activity: walking, running, driving, ...
- Sleep: rest/wake, in/out of bed, ...

Walking vs. time-of-day




Stage 2: Feature extraction

- Walking: cadence, stride-variability, asymmetry, ...
- Sleeping: time in bed, fragmentation, variability, ...



Stage 3: Feature Grouping



- Example: grouping of gait characteristics



Gait in Mild Alzheimer's Disease.
Feasibility of Multi-Center Measurement
in the Clinic and Home with Body-Worn Sensors.
A Pilot Study, Zetterberg et al, 2019, JAD



Sensor fusion

HYPNOS – Monitoring of type 2 diabetes patients



Macro-level: Diurnal Patterns

- ▶ Estimated functional mean, bmi, and age effects



Activity Pattern = mean pattern + Age*(Age Effect) + BMI*(BMI effect) + other effects